

Antibodies Product Datasheet

Product Name:	PPP2R3C Antibody
Catalog Number:	R12-3355
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	PPP2R3C Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP2R3C.
Concentration:	1 mg/ml
Gene Name:	PPP2R3C
Other Name:	Protein phosphatase subunit G5PR; Rhabdomyosarcoma antigen MU-RMS-40.6A; Rhabdomyosarcoma antigen Mu-RMS-40.6C
Human Swiss Prot No.:	Q969Q6
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	PPP2R3C Antibodyantibody detects endogenous levels of PPP2R3C.
Storage and Stability:	Stable at -20°C for at least 1 year.