

Antibodies Product Datasheet

Product Name:	RAR β Antibody
Catalog Number:	C0315
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Retinoic Acid Receptor beta. AA range:331-380
Concentration:	1 mg/ml
Gene Name:	RARB
Protein Name:	Retinoic acid receptor beta
Other Name:	RARB; HAP; NR1B2; Retinoic acid receptor beta; RAR-beta; HBV-activated protein; Nuclear receptor subfamily 1 group B member 2; RAR-epsilon
Human Swiss Prot No.:	P10826
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RAR β Polyclonal Antibody detects endogenous levels of RAR β protein.
Storage and Stability:	-20°C/1 year