

Antibodies Product Datasheet

Product Name:	ALS2 Antibody
Catalog Number:	C42550
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 390-470
Concentration:	1 mg/ml
Gene Name:	ALS2 ALS2CR6 KIAA1563
Protein Name:	Alsin (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 6 protein) (Amyotrophic lateral sclerosis 2 protein)
Other Name:	
Human Swiss Prot No.:	Q96Q42
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ALS2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year