

Antibodies Product Datasheet

Product Name:	Axl (Phospho-Tyr691) Antibody
Catalog Number:	A8079
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AXL around the phosphorylation site of Tyr691. AA range:657-706
Concentration:	1 mg/ml
Gene Name:	AXL
Protein Name:	Tyrosine-protein kinase receptor UFO
Other Name:	AXL; UFO; Tyrosine-protein kinase receptor UFO; AXL oncogene
Human Swiss Prot No.:	P30530
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Axl (Y691) Polyclonal Antibody detects endogenous levels of Axl protein only when phosphorylated at Y691.
Storage and Stability:	-20°C/1 year