

Antibodies Product Datasheet

Product Name:	UBL4A Antibody
Catalog Number:	C51159
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal UBL4A
Concentration:	1 mg/ml
Gene Name:	UBL4A DXS254E GDX UBL4
Protein Name:	Ubiquitin-like protein 4A (Ubiquitin-like protein GDX)
Other Name:	Ubiquitin-like protein 4A (Ubiquitin-like protein GDX)
Human Swiss Prot No.:	P11441
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UBL4A at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year