

Antibodies Product Datasheet

Product Name:	GABA A Receptor γ 2 Antibody
Catalog Number:	C40056
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GABA A Receptor γ 2
Concentration:	1 mg/ml
Gene Name:	GABRG2
Protein Name:	Gamma-aminobutyric acid receptor subunit gamma-2 (GABA(A) receptor subunit gamma-2)
Other Name:	Gamma-aminobutyric acid receptor subunit gamma-2 (GABA(A) receptor subunit gamma-2); YM3554
Human Swiss Prot No.:	P18507
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABA A Receptor γ 2 protein(A229) detects endogenous levels of GABA A Receptor γ 2
Storage and Stability:	-20°C/1 year