

Antibodies Product Datasheet

Product Name:	DDR1 Antibody
Catalog Number:	R12-2110
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DDR1. AA range:749-798
Concentration:	1 mg/ml
Gene Name:	DDR1
Protein Name:	Epithelial discoidin domain-containing receptor 1
Other Name:	DDR1; CAK; EDDR1; NEP; NTRK4; PTK3A; RTK6; TRKE; Epithelial discoidin domain-containing receptor 1; Epithelial discoidin domain receptor 1; CD167 antigen-like family member A; Cell adhesion kinase; Discoidin receptor tyrosine kinase; HGK2;
Human Swiss Prot No.:	Q08345
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DDR1 Polyclonal Antibody detects endogenous levels of DDR1 protein.
Storage and Stability:	-20°C/1 year