

Antibodies Product Datasheet

Product Name:	p120 (Phospho-Tyr228) Antibody
Catalog Number:	A0891
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Catenin-delta1 around the phosphorylation site of Tyr228. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	CTNND1
Protein Name:	Catenin delta-1
Other Name:	CTNND1; KIAA0384; Catenin delta-1; Cadherin-associated Src substrate; CAS; p120 catenin; p120(ctn); p120(cas)
Human Swiss Prot No.:	O60716
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-p120 (Y228) Polyclonal Antibody detects endogenous levels of p120 protein only when phosphorylated at Y228.
Storage and Stability:	-20°C/1 year