

Antibodies Product Datasheet

Product Name:	CYP21A2 Antibody
Catalog Number:	C12252
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 21A2. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	CYP21A2
Protein Name:	Steroid 21-hydroxylase
Other Name:	CYP21A2; CYP21; CYP21B; Steroid 21-hydroxylase; 21-OHase; Cytochrome P-450c21; Cytochrome P450 21; Cytochrome P450 XXI; Cytochrome P450-C21; Cytochrome P450-C21B
Human Swiss Prot No.:	P08686
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP21A2 Polyclonal Antibody detects endogenous levels of CYP21A2 protein.
Storage and Stability:	-20°C/1 year