

## Antibodies Product Datasheet

<b>Product Name:</b>	Cytochrome P450 24A1 Antibody
<b>Catalog Number:</b>	R12-2674
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Cytochrome P450 24A1 Antibody
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit Polyclonal
<b>Purification:</b>	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 24A1.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CYP24A1
<b>Other Name:</b>	1,25-dihydroxyvitamin D(3) 24-hydroxylase mitochondrial; Vitamin D(3) 24-hydroxylase; 24-Ohase; Cytochrome P450 24A1; P450-CC24; CYP24A1; CYP24
<b>Human Swiss Prot No.:</b>	Q07973
<b>Formulation:</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Specificity:</b>	Cytochrome P450 24A1 Antibodyantibody detects endogenous levels of Cytochrome P450 24A1.
<b>Storage and Stability:</b>	Stable at -20°C for at least 1 year.