

Antibodies Product Datasheet

Product Name:	PDZD3 Antibody
Catalog Number:	C50556
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PDZD3
Concentration:	1 mg/ml
Gene Name:	PDZD3 IKEPP NHERF4 PDZK2 DLNB27
Protein Name:	Na(+)/H(+) exchange regulatory cofactor NHE-RF4 (NHERF-4) (Intestinal and kidney-enriched PDZ protein) (Natrium-phosphate cotransporter IIa C-terminal-associated protein 2) (Na/Pi cotransporter C-terminal-associated protein 2) (NaPi-Cap2) (PDZ domain-cont
Other Name:	
Human Swiss Prot No.:	Q86UT5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PDZD3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year