

Antibodies Product Datasheet

Product Name:	P5CR3 Antibody
Catalog Number:	C51175
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal P5CR3
Concentration:	1 mg/ml
Gene Name:	PYCRL
Protein Name:	Pyrroline-5-carboxylate reductase 3 (P5C reductase 3) (P5CR 3) (EC 1.5.1.2) (Pyrroline-5-carboxylate reductase-like protein)
Other Name:	Pyrroline-5-carboxylate reductase 3 (P5C reductase 3) (P5CR 3) (EC 1.5.1.2) (Pyrroline-5-carboxylate reductase-like protein)
Human Swiss Prot No.:	Q53H96
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of P5CR3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year