

Antibodies Product Datasheet

Product Name:	UHRF1 Monoclonal Antibody
Catalog Number:	N0763
Reactivity:	Human
Applications:	WB,IF/ICC
Dilutions:	Western Blot: 1/1000 - 1/2000. Immunofluorescence: 1/100 - 1/500. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant human UHRF1 (N-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	UHRF1
Protein Name:	E3 ubiquitin-protein ligase UHRF1
Other Name:	UHRF1; ICBP90; NP95; RNF106; E3 ubiquitin-protein ligase UHRF1; Inverted CCAAT box-binding protein of 90 kDa; Nuclear protein 95; Nuclear zinc finger protein Np95; HuNp95; hNp95; RING finger protein 106; Transcription factor ICBP90; Ubiquit
Human Swiss Prot No.:	Q96T88
Formulation:	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Specificity:	UHRF1 Monoclonal Antibody detects endogenous levels of UHRF1 protein.
Storage and Stability:	-20°C/1 year