

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 10 (ABT056) Mouse Monoclonal Antibody
Catalog Number:	V0282
Clone ID:	ABT056
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT10
Gene Name:	P13645
Other Name:	BCIE, BIE, CK 10, CK-10, Cytokeratin-10, EHK, K10, K1C10_HUMAN, Keratin 10, Keratin 10 type I, Keratin, Keratin type I cytoskeletal 10, Keratin type I cytoskeletal 59 kDa, Keratin-10, Keratin10, KPP, KRT10, type I cytoskeletal 10
Isotype	IgG1, Kappa
Localization:	Skin
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK10 protein, and shows no cross reaction with CK4, 5, 6, 7, 8, 14, 15, 18, 19.
Storage and Stability:	-20°C/1 year