

Antibodies Product Datasheet

Product Name:	Anti-Human Uroplakin III (ABT523) Mouse Monoclonal Antibody
Catalog Number:	V0473
Clone ID:	ABT523
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	UPK3A
Gene Name:	O75631
Other Name:	MGC119178, OTTHUMP00000028951, UP III, UP3a, UPIII, UPIIIA, UPK 3, UPK 3A, UPK3, UPK3A, UPK3A_HUMAN, Uroplakin 3A, Uroplakin III, Uroplakin-3a, Uroplakin3A, UroplakinIII
Isotype	IgG1, kappa
Localization:	Bladder
Positive Tissue Control:	Membranous, Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Uroplakin III protein.
Storage and Stability:	-20°C/1 year