
Antibodies Product Datasheet

Product Name:	Claudin-2 Mouse Monoclonal Antibody(8H11)
Catalog Number:	N1722
Reactivity:	Human, Mouse, Rat
Applications:	WB,IHC
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Claudin-2
Concentration:	1 mg/ml
Gene Name:	CLDN2 PSEC0059 SP82 UNQ705/PRO1356
Protein Name:	Claudin-2
Other Name:	Claudin-2 (SP82)
Human Swiss Prot No.:	P57739
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Claudin-2
Storage and Stability:	-20°C/1 year