

Antibodies Product Datasheet

Product Name:	SGK1 Antibody
Catalog Number:	B8475
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from SGK1 . at AA range: 180-260
Concentration:	1 mg/ml
Gene Name:	SGK1
Protein Name:	Serine/threonine-protein kinase Sgk1
Other Name:	SGK1; SGK; Serine/threonine-protein kinase Sgk1; Serum/glucocorticoid-regulated kinase 1
Human Swiss Prot No.:	O00141
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SGK1 Polyclonal Antibody detects endogenous levels of SGK1 protein.
Storage and Stability:	-20°C/1 year