
Antibodies Product Datasheet

Product Name:	Aladin Antibody
Catalog Number:	C14368
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Aladin . at AA range: 360-440
Concentration:	1 mg/ml
Gene Name:	AAAS
Protein Name:	Aladin
Other Name:	AAAS; ADRACALA; GL003; Aladin; Adracalin
Human Swiss Prot No.:	Q9NRG9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Aladin Polyclonal Antibody detects endogenous levels of Aladin protein.
Storage and Stability:	-20°C/1 year