

Antibodies Product Datasheet

Product Name:	RBAP48 Antibody
Catalog Number:	C48395
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RBAP48
Concentration:	1 mg/ml
Gene Name:	RBBP4 RBAP48
Protein Name:	Histone-binding protein RBBP4 (Chromatin assembly factor 1 subunit C) (CAF-1 subunit C) (Chromatin assembly factor I p48 subunit) (CAF-I 48 kDa subunit) (CAF-I p48) (Nucleosome-remodeling factor subunit RBAP48) (Retinoblastoma-binding protein 4) (RBBP-4)
Other Name:	
Human Swiss Prot No.:	Q09028
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RBAP48 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year