

---

## Antibodies Product Datasheet

<b>Product Name:</b>	FAK Antibody
<b>Catalog Number:</b>	R12-2762
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC, ELISA
<b>Dilutions:</b>	FAK Antibody
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit Polyclonal
<b>Purification:</b>	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human FAK.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PTK2
<b>Other Name:</b>	FADK 1; FAK1; Focal adhesion kinase 1; pp125FAK; Protein- tyrosine kinase 2; PTK2
<b>Human Swiss Prot No.:</b>	Q05397
<b>Formulation:</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Specificity:</b>	FAK Antibodyantibody detects endogenous levels of FAK.
<b>Storage and Stability:</b>	Stable at -20°C for at least 1 year.