

Antibodies Product Datasheet

Product Name:	Flk-1/VEGFR2 Monoclonal Antibody
Catalog Number:	N0276
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant extracellular fragment of human Flk-1 (aa20-764) fused with hlgGFc tag expressed in HEK293 cells.
Concentration:	
Gene Name:	KDR
Protein Name:	Vascular endothelial growth factor receptor 2
Other Name:	KDR; FLK1; VEGFR2; Vascular endothelial growth factor receptor 2; VEGFR-2; Fetal liver kinase 1; FLK-1; Kinase insert domain receptor; KDR; Protein-tyrosine kinase receptor flk-1; CD antigen CD309
Human Swiss Prot No.:	P35968
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Flk-1 Monoclonal Antibody detects endogenous levels of Flk-1 protein.
Storage and Stability:	-20°C/1 year