

Antibodies Product Datasheet

Product Name:	PLEKHB1 Antibody
Catalog Number:	C50048
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PLEKHB1
Concentration:	1 mg/ml
Gene Name:	PLEKHB1 EVT1 KPL1 PHR1 PHRET1
Protein Name:	Pleckstrin homology domain-containing family B member 1 (PH domain-containing family B member 1) (Evectin-1) (PH domain-containing protein in retina 1) (PHRET1) (Pleckstrin homology domain retinal protein 1)
Other Name:	
Human Swiss Prot No.:	Q9UF11
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PLEKHB1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year