

Antibodies Product Datasheet

Product Name:	Cleaved-Ephrin-A4 (S170) Antibody
Catalog Number:	L0397
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Cleaved-Ephrin-A4 (S170)
Concentration:	1 mg/ml
Gene Name:	EFNA4 EPLG4 LERK4
Protein Name:	Cleaved-Ephrin-A4 (S170)
Other Name:	Ephrin-A4 (EPH-related receptor tyrosine kinase ligand 4) (LERK-4)
Human Swiss Prot No.:	P52798
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse Cleaved-Ephrin-A4 (S170, protein was cleaved amino acid sequence between 170-171)
Storage and Stability:	-20°C/1 year