

Antibodies Product Datasheet

Product Name:	MDMX Antibody
Catalog Number:	B8369
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MDM4. AA range:336-385
Concentration:	1 mg/ml
Gene Name:	MDM4
Protein Name:	Protein Mdm4
Other Name:	MDM4; MDMX; Protein Mdm4; Double minute 4 protein; Mdm2-like p53-binding protein; Protein Mdmx; p53-binding protein Mdm4
Human Swiss Prot No.:	O15151
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MDMX Polyclonal Antibody detects endogenous levels of MDMX protein.
Storage and Stability:	-20°C/1 year