
Antibodies Product Datasheet

Product Name:	ASK1 (Phospho-Thr838) Antibody
Catalog Number:	A9432
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ASK1 (Phospho Thr838)
Concentration:	1 mg/ml
Gene Name:	MAP3K5 ASK1 MAPKKK5 MEKK5
Protein Name:	ASK1 (Phospho Thr838)
Other Name:	Mitogen-activated protein kinase kinase kinase 5 (EC 2.7.11.25) (Apoptosis signal-regulating kinase 1) (ASK-1) (MAPK/ERK kinase kinase 5) (MEK kinase 5) (MEKK 5)
Human Swiss Prot No.:	Q99683
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human ASK1 (Phospho Thr838)
Storage and Stability:	-20°C/1 year