

Antibodies Product Datasheet

Product Name:	NUDC (Phospho-Ser326) Antibody
Catalog Number:	A8388
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NudC around the phosphorylation site of Ser326. AA range:282-331
Concentration:	1 mg/ml
Gene Name:	NUDC
Protein Name:	Nuclear migration protein nudC
Other Name:	NUDC; Nuclear migration protein nudC; Nuclear distribution protein C homolog
Human Swiss Prot No.:	Q9Y266
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-NUDC (S326) Polyclonal Antibody detects endogenous levels of NUDC protein only when phosphorylated at S326.
Storage and Stability:	-20°C/1 year