

Antibodies Product Datasheet

Product Name:	Emx2 Antibody
Catalog Number:	C10877
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EMX2. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	EMX2
Protein Name:	Homeobox protein EMX2
Other Name:	EMX2; Homeobox protein EMX2; Empty spiracles homolog 2; Empty spiracles-like protein 2
Human Swiss Prot No.:	Q04743
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Emx2 Polyclonal Antibody detects endogenous levels of Emx2 protein.
Storage and Stability:	-20°C/1 year