
Antibodies Product Datasheet

Product Name:	NGF Antibody
Catalog Number:	C30199
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human NGF. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	NGF
Protein Name:	Beta-nerve growth factor
Other Name:	NGF; NGFB; Beta-nerve growth factor; Beta-NGF
Human Swiss Prot No.:	P01138
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NGF Polyclonal Antibody detects endogenous levels of NGF protein.
Storage and Stability:	-20°C/1 year