

## Antibodies Product Datasheet

<b>Product Name:</b>	MST3 (Phospho-Thr18) Antibody
<b>Catalog Number:</b>	A8150
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human MST3 (Phospho-Thr18)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	STK24 MST3 STK3
<b>Protein Name:</b>	MST3 (Phospho-Thr18)
<b>Other Name:</b>	Serine/threonine-protein kinase 24 (EC 2.7.11.1) (Mammalian STE20-like protein kinase 3) (MST-3) (STE20-like kinase MST3) [Cleaved into: Serine/threonine-protein kinase 24 36 kDa subunit (Mammalian STE20-like protein kinase 3 N-terminal) (MST3/N); Serine/
<b>Human Swiss Prot No.:</b>	Q9Y6E0
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.203% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of MST3 (Phospho-Thr18) at Human, Mouse, Rat
<b>Storage and Stability:</b>	-20°C/1 year