

Antibodies Product Datasheet

Product Name:	ApoC-III Antibody
Catalog Number:	C30070
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human APOC3. AA range:46-95
Concentration:	1 mg/ml
Gene Name:	APOC3
Protein Name:	Apolipoprotein C-III
Other Name:	APOC3; Apolipoprotein C-III; Apo-CIII; ApoC-III; Apolipoprotein C3
Human Swiss Prot No.:	P02656
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ApoC-III Polyclonal Antibody detects endogenous levels of ApoC-III protein.
Storage and Stability:	-20°C/1 year