

## **Assay Biotechnology Company**

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

## **Antibodies Product Datasheet**

Product Name: p39 Antibody

Catalog Number: C10968

Reactivity: Human, Rat, Mouse

**Applications:** WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000. Not yet tested in other applications.

**Modification:** Unmodified/Total

**Source:** Polyclonal Rabbit

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Immunogen:** The antiserum was produced against synthesized peptide derived from

human CDK5R2. AA range:81-130

Concentration: 1 mg/ml

Gene Name: CDK5R2

**Protein Name:** Cyclin-dependent kinase 5 activator 2

Other Name: CDK5R2; NCK5AI; Cyclin-dependent kinase 5 activator 2; CDK5

activator 2; Cyclin-dependent kinase 5 regulatory subunit 2; p39; p391

**Human Swiss Prot No.:** Q13319

**Formulation:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

**Specificity:** p39 Polyclonal Antibody detects endogenous levels of p39 protein.

Storage and Stability: -20°C/1 year