

## Antibodies Product Datasheet

<b>Product Name:</b>	JAK3 (Phospho-Tyr981) Antibody
<b>Catalog Number:</b>	A9448
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human JAK3 (Phospho Tyr981)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	JAK3
<b>Protein Name:</b>	JAK3 (Phospho Tyr981)
<b>Other Name:</b>	Tyrosine-protein kinase JAK3 (EC 2.7.10.2) (Janus kinase 3) (JAK-3) (Leukocyte janus kinase) (L-JAK)
<b>Human Swiss Prot No.:</b>	P52333
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human, Mouse, Rat JAK3 (Phospho Tyr981)
<b>Storage and Stability:</b>	-20°C/1 year