

Antibodies Product Datasheet

Product Name:	WASP (Phospho-Tyr290) Antibody
Catalog Number:	A0597
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human WASP around the phosphorylation site of Tyr290. AA range:256-305
Concentration:	1 mg/ml
Gene Name:	WAS
Protein Name:	Wiskott-Aldrich syndrome protein
Other Name:	WAS; IMD2; Wiskott-Aldrich syndrome protein; WASp
Human Swiss Prot No.:	P42768
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-WASP (Y290) Polyclonal Antibody detects endogenous levels of WASP protein only when phosphorylated at Y290.
Storage and Stability:	-20°C/1 year