

Antibodies Product Datasheet

Product Name:	CCDC59 Antibody
Catalog Number:	C49926
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CCDC59
Concentration:	1 mg/ml
Gene Name:	CCDC59 BR22 TAP26 HSPC128
Protein Name:	Thyroid transcription factor 1-associated protein 26 (TTF-1-associated protein 26) (Coiled-coil domain-containing protein 59) (TTF-1-associated protein BR2)
Other Name:	
Human Swiss Prot No.:	Q9P031
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CCDC59 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year