

Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Rab 3 GAP p130 Antibody |
| Catalog Number: | C18124 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | WB, IHC-p, IF(paraffin section), ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human RAB3GAP1. AA range:538-587 |
| Concentration: | 1 mg/ml |
| Gene Name: | RAB3GAP1 |
| Protein Name: | Rab3 GTPase-activating protein catalytic subunit |
| Other Name: | RAB3GAP1; KIAA0066; RAB3GAP; Rab3 GTPase-activating protein catalytic subunit; RAB3 GTPase-activating protein 130 kDa subunit; Rab3-GAP p130; Rab3-GAP |
| Human Swiss Prot No.: | Q15042 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Rab 3 GAP p130 Polyclonal Antibody detects endogenous levels of Rab 3 GAP p130 protein. |
| Storage and Stability: | -20°C/1 year |