
Antibodies Product Datasheet

Product Name:	Dkk-3 Monoclonal Antibody
Catalog Number:	N0198
Reactivity:	Human, Monkey
Applications:	WB, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Dkk-3 expressed in E. Coli.
Concentration:	
Gene Name:	DKK3
Protein Name:	Dickkopf-related protein 3
Other Name:	DKK3; REIC; Dickkopf-related protein 3; Dickkopf-3; Dkk-3; hDkk-3
Human Swiss Prot No.:	Q9UBP4
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Dkk-3 Monoclonal Antibody detects endogenous levels of Dkk-3 protein.
Storage and Stability:	-20°C/1 year