
Antibodies Product Datasheet

Product Name:	eIF4B (Phospho-Ser406) Antibody
Catalog Number:	A9640
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human eIF4B (Ser406)
Concentration:	1 mg/ml
Gene Name:	EIF4B
Protein Name:	eIF4B (Ser406)
Other Name:	Eukaryotic translation initiation factor 4B (eIF-4B)
Human Swiss Prot No.:	P23588
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat eIF4B (phospho-Ser406)
Storage and Stability:	-20°C/1 year