

Antibodies Product Datasheet

Product Name:	CKR-4 Antibody
Catalog Number:	C12128
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCR4. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	CCR4
Protein Name:	C-C chemokine receptor type 4
Other Name:	CCR4; CMKBR4; C-C chemokine receptor type 4; C-C CKR-4; CC-CKR-4; CCR-4; CCR4; K5-5; CD antigen CD194
Human Swiss Prot No.:	P51679
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CKR-4 Polyclonal Antibody detects endogenous levels of CKR-4 protein.
Storage and Stability:	-20°C/1 year