Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: E-cadherin Antibody

Catalog Number: C0147

Reactivity: Human, Rat, Mouse

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in

other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human E-cadherin. AA range:1-50

Concentration: 1 mg/ml

Gene Name: CDH1

Protein Name: Cadherin-1

Other Name: CDH1; CDHE; UVO; Cadherin-1; CAM 120/80; Epithelial cadherin; E-

cadherin; Uvomorulin; CD antigen CD324

Human Swiss Prot No.: P12830

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: E-cadherin Polyclonal Antibody detects endogenous levels of E-cadherin

protein.

Storage and Stability: -20°C/1 year