

Antibodies Product Datasheet

Product Name:	Cytochrome P450 2C19 Antibody
Catalog Number:	R12-2675
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	Cytochrome P450 2C19 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 2C19.
Concentration:	1 mg/ml
Gene Name:	CYP2C19
Other Name:	Cytochrome P450 2C19; (R)-limonene 6-monooxygenase; (S)-limonene 6-monooxygenase; (S)-limonene 7-monooxygenase; CYP11C19; P450-11A; Mephenytoin 4-hydroxylase; CYP11C17; P450-254C; CYP2C19; CP2CJ
Human Swiss Prot No.:	P33261
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	Cytochrome P450 2C19 Antibody antibody detects endogenous levels of Cytochrome P450 2C19.
Storage and Stability:	Stable at -20°C for at least 1 year.