

Antibodies Product Datasheet

Product Name:	PBEF Monoclonal Antibody
Catalog Number:	N0510
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of PBEF (aa338-479) expressed in E. Coli.
Concentration:	
Gene Name:	NAMPT
Protein Name:	Nicotinamide phosphoribosyltransferase
Other Name:	NAMPT; PBEF; PBEF1; Nicotinamide phosphoribosyltransferase; NAmPRTase; Nampt; Pre-B-cell colony-enhancing factor 1; Pre-B cell-enhancing factor; Visfatin
Human Swiss Prot No.:	P43490
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	PBEF Monoclonal Antibody detects endogenous levels of PBEF protein.
Storage and Stability:	-20°C/1 year