

## Antibodies Product Datasheet

<b>Product Name:</b>	BRF1 Antibody
<b>Catalog Number:</b>	R12-2546
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IHC, ELISA
<b>Dilutions:</b>	BRF1 Antibody
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Rabbit Polyclonal
<b>Purification:</b>	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human BRF1.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	BRF1
<b>Other Name:</b>	Transcription factor IIIB 90 kDa subunit; TFIIIB90; hTFIIIB90; B-related factor 1; BRF-1; hBRF; TAF3B2; TATA box-binding protein-associated factor RNA polymerase III subunit 2; BRF1; BRF; GTF3B; TAF3B2; TAF3C
<b>Human Swiss Prot No.:</b>	Q92994
<b>Formulation:</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Specificity:</b>	BRF1 Antibody antibody detects endogenous levels of BRF1.
<b>Storage and Stability:</b>	Stable at -20°C for at least 1 year.