

Antibodies Product Datasheet

Product Name:	RIPK5 Antibody
Catalog Number:	C48928
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RIPK5
Concentration:	1 mg/ml
Gene Name:	DSTYK KIAA0472 RIP5 RIPK5 SGK496 HDCMD38P
Protein Name:	Dual serine/threonine and tyrosine protein kinase (EC 2.7.12.1) (Dusty protein kinase) (Dusty PK) (RIP-homologous kinase) (Receptor-interacting serine/threonine-protein kinase 5) (Sugen kinase 496) (SgK496)
Other Name:	
Human Swiss Prot No.:	Q6XUX3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RIPK5 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year