

Antibodies Product Datasheet

Product Name:	NDUFB4 Antibody
Catalog Number:	C50085
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NDUFB4
Concentration:	1 mg/ml
Gene Name:	NDUFB4
Protein Name:	NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 4 (Complex I-B15) (CI-B15) (NADH-ubiquinone oxidoreductase B15 subunit)
Other Name:	
Human Swiss Prot No.:	O95168
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NDUFB4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year