

Antibodies Product Datasheet

Product Name:	Dab1 (Phospho-Tyr232) Antibody
Catalog Number:	A7055
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Dab1 around the phosphorylation site of Tyr232. AA range:199-248
Concentration:	1 mg/ml
Gene Name:	DAB1
Protein Name:	Disabled homolog 1
Other Name:	DAB1; Disabled homolog 1
Human Swiss Prot No.:	O75553
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-Dab1 (Y232) Polyclonal Antibody detects endogenous levels of Dab1 protein only when phosphorylated at Y232.
Storage and Stability:	-20°C/1 year