

Antibodies Product Datasheet

Product Name:	SSTR2 Antibody
Catalog Number:	C30261
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SSTR2. AA range:220-270
Concentration:	1 mg/ml
Gene Name:	SSTR2
Protein Name:	somatostatin receptor 2
Other Name:	Somatostatin receptor type 2 (SS-2-R) (SS2-R) (SS2R) (SRIF-1)
Human Swiss Prot No.:	P30874
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SSTR2 Polyclonal Antibody detects endogenous levels of somatostatin receptor 2
Storage and Stability:	-20°C/1 year