

Antibodies Product Datasheet

Product Name:	ZC3H13 Antibody
Catalog Number:	R12-3687
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	ZC3H13 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZC3H13.
Concentration:	1 mg/ml
Gene Name:	ZC3H13
Other Name:	Zinc finger CCCH domain-containing protein 13; ZC3H13; KIAA0853
Human Swiss Prot No.:	Q5T200
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ZC3H13 Antibodyantibody detects endogenous levels of ZC3H13.
Storage and Stability:	Stable at -20°C for at least 1 year.