

Antibodies Product Datasheet

Product Name:	NY-CO-1 Antibody
Catalog Number:	C10140
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SDCG1. AA range:881-930
Concentration:	1 mg/ml
Gene Name:	NEMF
Protein Name:	Nuclear export mediator factor NEMF
Other Name:	NEMF; SDCCAG1; Nuclear export mediator factor NEMF; Antigen NY-CO-1; Serologically defined colon cancer antigen 1
Human Swiss Prot No.:	O60524
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NY-CO-1 Polyclonal Antibody detects endogenous levels of NY-CO-1 protein.
Storage and Stability:	-20°C/1 year