

## Antibodies Product Datasheet

<b>Product Name:</b>	AF9 Antibody
<b>Catalog Number:</b>	C42515
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human AF9
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MLLT3 AF9 YEATS3
<b>Protein Name:</b>	AF9
<b>Other Name:</b>	Protein AF-9 (ALL1-fused gene from chromosome 9 protein) (Myeloid/lymphoid or mixed-lineage leukemia translocated to chromosome 3 protein) (YEATS domain-containing protein 3)
<b>Human Swiss Prot No.:</b>	P42568
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.153% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of AF9 at Human/Mouse
<b>Storage and Stability:</b>	-20°C/1 year