

Antibodies Product Datasheet

Product Name:	GFP Monoclonal Antibody(Mix)
Catalog Number:	N1389
Reactivity:	Species independent
Applications:	WB,IF/ICC,IP
Dilutions:	WB: 1:5000 IP: 1:200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of GFP
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detectsThe antibody detects RFP, YFP, CFP.
Storage and Stability:	-20°C/1 year